INFORMATION INSERT

Ranfac Bone Marrow Biopsy Needle

Sterile: Contents sterile unless package is opened, damaged or expired.

A. INDICATIONS FOR USE

The Ranfac Bone Marrow Biopsy Needle is intended for the purpose of harvesting bone marrow specimens.

B. CONTRAINDICATIONS

For use only for biopsies of bone marrow as determined by a licensed physician. The device is intended to be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone marrow biopsy. The procedure should be performed on patients that are suitable for such procedure only.

WARNING: This device must not be used for sternal procedures.

CAUTION: For single patient use only. The design of this device may not perform as intended by the manufacturer if it is re-used. The manufacturer cannot guarantee the performance, safety and reliability of a reprocessed device.

After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose of in accordance with applicable laws and regulations. Carefully place the used needle in a sharps biohazard container after the procedure is completed.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

C. DIRECTIONS FOR USE

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the contents from the tray/pouch.
3. Following proper placing technique, locate the needle (1, 2) adjacent to the biopsy site.
4. Using gentle, but firm, pressure advance the needle (1, 2) by rotating it in an alternating clockwise - counterclockwise motion. Entrance into the marrow cavity is generally detected by decreased resistance.
5. Unlock the stylet cap (2) by rotating it counterclockwise, and remove the stylet (2) from the cannula (1).
6. Use the needle vise (5) provided to cap the stylet (2) after use for safe disposal.
7. Very slowly and gently, advance the needle (1) using clockwise - counterclockwise motion until an adequate depth is achieved (2-3cm).
8. Insert the marked probe (3) without force to check sample length in needle lumen. Sample length is referenced by the marks in comparison to the needle hub. Remove the probe.
9. Deploy the gray lever by moving it from the "OPEN" to the "CLOSED" position.
10. Slowly withdraw the needle (1) from the patient.
11. Rotate the gray lever from the "CLOSED" to the "OPEN" position.
12. Place the probe guide (4) onto the distal tip of the cannula (1) as shown. Insert the probe (3) into the distal end of the probe guide (4). Advance the probe (3) into the needle (1) and collect the biopsy specimen as it emerges from the opening in the handle (proximal end).
13. Use the needle vise (5) provided to cap the needle (1) after use for safe disposal.


CAUTION: If any difficulty is encountered during the transfer of the sample, DO NOT force the sample. Remove the probe (3), activate the grey lever to the "CLOSED" position and return it to the "OPEN" position. Re-insert the probe (3) and transfer the specimen out the proximal end of the needle (1)

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

U.S. Pats. 5,522,398 and 5,634,473

1-800-2RANFAC (272-6322)

Ranfac Corp.
30 Doherty Ave.
Avon, MA 02322-0635
Tel: (800) 2-RANFAC
(508) 588-4400
Fax: (508) 584-8588
Email: info@ranfac.com
www.ranfac.com

Distributed By:

CoreLink, LLC
888.349.7808
CoreLinkSurgical.com

P/N 09653-01P REV C

Bone Marrow Aspiration Needle

Sterile: Contents sterile unless package is opened, damaged or expired.

A. INDICATIONS FOR USE

The Bone Marrow Aspiration Needle is intended for use in aspirating a bone marrow specimen

B. CONTRAINDICATIONS

Use only for bone marrow aspiration as determined by a licensed physician. The device is intended to be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone marrow aspiration. The procedure should be performed on patients that are suitable for such procedure only.

CAUTION: For single patient use only. The design of this device may not perform as intended by the manufacturer if it is re-used. The manufacturer cannot guarantee the performance, safety and reliability of a reprocessed device.

After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose of in accordance with applicable laws and regulations. Carefully place the used needle in a sharps biohazard container after the procedure is completed.

Note: *These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.*

C. DIRECTIONS FOR USE

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the needle from the tray/pouch and inspect the cutting edge of both cannula and stylet for any damage or other imperfections that would prevent the proper operation of the needle.

Note: Confirm that the match ground points of the stylet and needle are aligned properly by matching the pad print markings (if present) on the same side of the needle.

3. Following proper placing technique, locate the needle adjacent to the aspiration site.
4. Using gentle, but firm pressure, advance the needle by rotating in an alternating clockwise-counterclockwise motion. Entrance to the marrow cavity is generally detected by decreased resistance.
5. Rotate the stylet cap counterclockwise and pull the stylet out of the cannula.
6. Attach a luer slip syringe to the cannula hub

Note: *Avoid using glass syringes because they may result in an inadequate seal and may also break or fragment.*

7. Apply negative pressure by quickly withdrawing the syringe plunger. Disengage the syringe and remove the aspirated sample.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

1-800-2RANFAC (272-6322)

Ranfac Corp.
30 Doherty Ave.
Aven. MA 02322-0635
Tel: (800) 2-RANFAC
(508) 588-4400
Fax: (508) 584-8588
Email: info@ranfac.com
www.ranfac.com

Distributed By:

CoreLink, LLC
888.349.7808
CoreLinkSurgical.com

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ENGLISH / US

Stérile: Contents sterile unless package is opened, damaged or expired.

A. INDICATIONS FOR USE

The Ranfac Snarecoil® Bone Marrow Biopsy Needle is intended for the purpose of harvesting bone marrow specimens.

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C. DIRECTIONS FOR USE

1. Inspect the package for proper integrity. If undamaged, open the package using aseptic technique.
2. Remove the contents from the tray/pouch.
3. Following proper placing technique, locate the needle adjacent to the biopsy site.
4. Using gentle, but firm, pressure advance the needle by rotating it in an alternating clockwise - counterclockwise motion. Entrance into the marrow cavity is generally detected by decreased resistance.
5. Unlock the stylet cap by rotating it counterclockwise, and remove the stylet from the cannula.
6. Use the needle vise provided to cap the stylet after use for safe disposal.
7. Very slowly and gently, advance the needle using clockwise - counterclockwise motion until an adequate depth is achieved (2-3cm).
8. Insert the marked probe without force to check sample length in needle lumen. Sample length is referenced by the marks in comparison to the needle hub. Remove the probe.
9. Deploy the gray lever by moving it from the "OPEN" to the "CLOSED" position.
10. Slowly withdraw the needle from the patient.
11. Rotate the gray lever from the "CLOSED" to the "OPEN" position.
12. Place the probe guide onto the distal tip of the cannula. Insert the probe into the distal end of the probe guide. Advance the probe into the needle and collect the biopsy specimen as it emerges from the opening in the handle (proximal end).
13. Use the needle vise provided to cap the needle after use for safe disposal.

CAUTION: If any difficulty is encountered during the transfer of the sample, DO NOT force the sample. Remove the probe, activate the gray lever to the "CLOSED" position and return it to the "OPEN" position. Re-insert the probe and transfer the specimen out the proximal end of the needle.

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